



# AVANTBIO™

***Chemically-defined, Animal Origin-free Cell Culture  
For Regenerative Medicine and Other Applications***



**Historically, Many Cell- and Tissue-based Therapy Protocols Utilize Cell Culture Media that Includes the Use of Components that are Originated from Human or Non-human Animal (e.g. bovine, porcine, rodent) Sources. These Animal-derived Components Include:**

- *Blood-derived Serum and Plasma*
- *Serum-derived Albumin, Transferrin and Fetuin Protein*
- *Replication Incompetent Feeder Layer Cells (both mortal and immortal cell lines)*
- *Pituitary Gland Tissue Extract*
- *Placental Amniotic Membranes*
- *Extracellular Matrix Proteins*
- *Blood-derived Platelet Lysate*

## Motivation for Cell- and Engineered Tissue-based Therapy Companies to Adopt Chemically-defined Animal Origin-free (cdAOF) Cell Culture Systems:

- *Lessens the risk of animal component originated adverse events, that include the transfer of animal-originated pathogens (e.g., mad cow/BSE, viruses) to humans as well as inflammatory (allergic) reactions to animal-derived cell culture media components*
- *Reduces the potential for cell and tissue therapy product variability by using quality-tested, controllable, chemically defined cell culture components (i.e. precise control of cell and tissue culture manufacturing conditions)*
- *Pressure from regulatory authorities (FDA, EMA) to pursue cdAOF Cell Culture Options*
- *Can reduce the cost of manufacturing therapeutic cells and tissue, especially in light of recent qualified fetal bovine serum (FBS) shortages as well as the increasing cost for qualified FBS and other animal components*
- *Pressure from animal welfare groups and the general public to employ cdAOF cell culture technology, as it is also “Animal Cruelty-free”*
- *Provides a competitive marketing advantage that highlights complete “animal cruelty-free” (e.g. fetal bovine serum-free) cell and tissue culture environment when carrying out the in vitro testing and marketing of OTC cosmeceuticals and other consumer personal care products*

**To Avoid the Use of Animal Derived Components for The Culture of Human Skin-derived Dermal Fibroblasts, Avantbio Developed an Improved, Chemically-defined Animal Origin-free (cdAOF) Cell Culture Supplement System That Does Not Include any Human or Other Animal-derived Products. The Following is a Collection of Representative Research and Development Results Derived from the Testing of This New Clinically-relevant cdAOF Cell Culture System.**



**HFSdaFREE KIT**  
**3-component supplement kit:**  
**(HFSdaFREE, HFGE &HFGE2)**



## ***Chemically-defined, Animal Origin-free (cdAOF) Cell and Tissue Culture Systems***

**Slides 6-8** Serial Propagation of Serum-starved Neonatal Human Dermal Fibroblasts (HDFn)

**Slides 9-10** Serial Propagation of cdAOF Primary Adult Human Dermal Fibroblasts (HDFa)


**Slides 11-12** Performance Test Establishing the Stability of the HFSdaFREE KIT Supplement System during – 20 °C Storage

**Slides 13-15** Propagation of HDFn on a Novel Bio-compatible Substratum: Silk-derived Fibroin Discs (Mats)

**Slides 16-18** Serial Propagation of a Serum-starved Normal Diploid Human Fetal Lung Fibroblast Line (WI-38)

**Slides 19-21** Serial Propagation of Serum-starved Human Hair Follicle Dermal Papilla Cells (HFDPC)


**Slides 22-23** Serial Propagation of Serum-starved Adipose-derived Human Mesenchymal Stem Cells (HMSC)



# Human Dermal Fibroblasts (HDF)

## *Applications in Basic and Preclinical Laboratory Research*

- ***2D Human Fibroblast and Other Mesenchymal Cell Cultures***
  - *Skin-derived Dermal fibroblasts*
  - *Cornea-derived Stromal Fibroblasts (Keratocytes)*
  - *Dermal Fibroblasts as Precursor Cells for Induced Pluripotent Stem Cells (iPSCs)*
- ***3D Bioprinted, Automated or Manually Reconstructed Connective Tissue, Using Human Fibroblasts***
  - *Connective Tissue (e.g., Skin Dermis, Corneal Stroma)*
  - *Other Connective Tissues*
  - *Tissue on a Chip (e.g., Skin Dermis, Corneal Stroma)*






# For Human Dermal Fibroblasts (HDF)

## *Clinical Applications in Regenerative and Esthetic Medicine Using Dermal Fibroblasts and Other Mesenchymal Cells*


- **2D Cell-based Therapies (e.g. Cutaneous, Cornea, Other)**
- **2D Cell-derived Therapeutic Products (e.g., Conditioned Medium and Exosomes)**
- **3D Reconstructed Connective Tissue Therapies (e.g. Cutaneous, Cornea, Other)**

## *Cosmeceutical Applications Using Dermal Fibroblasts and Other Mesenchymal Cells*

- **Cell-derived Cosmeceutical Products (e.g., Conditioned Medium and Exosomes)**



# The HFSdaFREE KIT Supplement System Supports the Efficient Post-primary Serial Propagation of Serum-starved Neonatal Human Dermal Fibroblasts (HDFn) Under Chemically-defined Animal Origin-free (cdAOF) Conditions




**Methods:** Cryopreserved post-primary serum-starved passage 1 neonatal human dermal fibroblast (HDFn) reared in medium m106 plus LSGS (Life Technologies; 2% FBS + growth factors), at the primary culture level, were placed into passage 2 cell culture, under chemically defined animal origin-free (cdAOF) cell culture conditions. Briefly, using AvantBio's HFSdaFREE KIT supplements in a 50/50 blend of EpiLife basal medium/Medium m106, cdAOF HDFn were serially passaged through passage 5 and then cryopreserved in cdAOF CRYOVIVE 5% DMSO cryofluid. The passage 5 HDFn were later thawed and plated into passage 6 culture at 2,500 cells/cm<sup>2</sup>, into rh Collagen-1 precoated 6 well plates, in triplicate, using HFSdaFREE KIT supplements as described above, or control HDFn cell cultures propagated on untreated cell wells, using a conventional animal product containing medium (EMEM + 10% FBS). HDFn's were serially passaged through passage 11, in the cdAOF and control animal product-containing culture conditions, as described above. Mean cumulative population doublings were calculated at the end of each passage, from triplicate wells, +/- the SEM. LSGS, EpiLife and medium m106 are trademarks of Life Technologies Corporation.




# Efficient Recovery and Serial Propagation of Post-primary cdAOF Adult\* Human Dermal Fibroblasts (HDFa), After 16 months in Cryopreservation. HDFa Were Placed Into Post-Primary Passage 1 Culture, Under Chemically-defined Animal Origin-free (cdAOF) Cell Culture Conditions, Using the HFSdaFREE KIT Supplement System (\*53 year old female)

Post-thaw Viability of Primary HDFa Cultured Under cdAOF Conditions




Serial Propagation of Primary HDFa Under cdAOF Conditions



**Methods:** Chemically defined animal origin-free (cdAOF) adult human dermal fibroblast (HDFa) primary cell cultures were established using HFSdaFREE KIT supplements and 50/50 blend of EpiLife basal medium and Medium m106, dissociated and then cryopreserved in cdAOF CRYOVIVE 5% DMSO cryofluid. After 16 months in liquid nitrogen vapor storage, the HDFa were plated into passage 1 culture at 2,500 cells/cm<sup>2</sup>, into rh Collagen-1 precoated 6 well plates, in triplicate. Passage 1 cells were serially passaged in the cdAOF culture environment using AvantBio's HFSdaFREE KIT supplements and the basal medium described above. Mean cumulative population doublings were calculated at the end of each passage, from triplicate wells, +/- the SEM. EpiLife and medium m106 are trademarks of Life Technologies Corporation.


**Efficient Recovery and Serial Propagation of Post-primary cdAOF Adult\* Human Dermal Fibroblasts (HDFa), After 16 months in Cryopreservation. HDFa Were Placed Into Post-Primary Passage 1 Culture, Under Chemically-defined Animal Origin-free (cdAOF) Cell Culture Conditions, Using the HFSdaFREE KIT Supplement System (\*53 year old female)**




**Passage 1, Day 1**



**Passage 1, Day 3**



**Passage 1, Day 5**




high magnification

**Passage 1, Day 5**


**Methods:** Chemically defined animal origin-free (cdAOF) adult human dermal fibroblast (HDFa) primary cell cultures were established from punch biopsies using HFSdaFREE KIT supplements (HFSdaFREE, HFGE, HFGE2) and 50/50 blend of EpiLife basal medium and Medium m106, dissociated and then cryopreserved in cdAOF CRYOVIVE 5% DMSO cryofluid. After 16 months in liquid nitrogen vapor storage, the HDFa were plated into passage 1 culture at 2,500 cells/cm<sup>2</sup>, into rh Collagen-1 precoated 6 well plates, in triplicate. Passage 1 cells were serially passaged in the cdAOF culture environment using AvantBio's HFSdaFREE KIT supplements and the basal medium described above. Representative photomicrographs were taken of cdAOF HDFa during passage 1 at days 1,3 and 5. EpiLife and medium m106 are trademarks of Life Technologies Corporation.

# Stability Performance Test on Neonatal Dermal Fibroblasts (HDFn): The HFSdaFREE Kit Supplement System is Stable When Stored at -20 °C, for Up to 3.6 years


Day 8: Mean Cell Densities




Day 8: Fixed & Stained HDFn



HFSdaFREE lot 76 + HFGE lot 77 + HFGE2 lot 78: 1.7 years in storage



HFSdaFREE lot 53 + HFGE lot 54 + HFGE2 lot 55: 2.2 years in storage



HFSdaFREE lot 33 + HFGE lot 34 + HFGE2 lot 35 : 3.6 years in storage


**Methods:** As described for slide 4, serum-starved, post-primary passage 2 human neonatal dermal fibroblast (HDFn) were cultured under chemically defined animal origin-free (cdAOF) conditions using the indicated lots of HFSdaFREE KIT supplements. HDFn were cultured on human recombinant Collagen-1 pre-coated cell wells (6 well format). Results are displayed as final mean cell densities after 8 days of culture, +/- the SEM. Seeding density = 250 cells/cm<sup>2</sup> (see red arrow). Representative stained cell images were taken at day 8. Pop. Dbl./day = average population doublings per day (a measure of relative growth rate). EpiLife Medium (ELM) and medium m106 are registered trademarks of Life Technologies Corporation.




# Stability Performance Test on Neonatal Dermal Fibroblasts (HDFn): The HFSdaFREE Kit Supplement System is Stable When Stored at -20 °C, for Up to 3.6 years

## Day 1 vs Day 8 Photomicrographs





HFSdaFREE lot 76 + HFGE lot 77 +  
HFGE2 lot 78: 1.7 years in storage



HFSdaFREE lot 53 + HFGE lot 54 +  
HFGE2 lot 55: 2.2 years in storage



HFSdaFREE lot 33 + HFGE lot 34 +  
HFGE2 lot 35 : 3.6 years in storage



**Methods:** As described for slide 4, serum-starved, post-primary passage 2 human neonatal dermal fibroblast (HDFn) were cultured under chemically defined animal origin-free (cdAOF) conditions using the indicated lots of HFSdaFREE KIT supplements. HDFn were cultured on human recombinant Collagen-1 pre-coated cell wells (6 well format). Results are displayed as final mean cell densities after 8 days of culture, +/- the SEM. Seeding density = 250 cells/cm<sup>2</sup> (see red arrow). Representative stained cell images were taken at day 8. EpiLife Medium (ELM) and medium m106 are registered trademarks of Life Technologies Corporation.

# Bio-compatible Engineered Silk-derived Fibroin Protein Can Be Fabricated Into a Variety Configurations, in Some cases for Applications in Regenerative Medicine, That Include the Delivery of Therapeutic Cells or Tissue on Discs (Mats), Scaffolds (Fibers) and Hydrogels

International Journal of  
Molecular Sciences

Review

## Silk Fibroin as a Functional Biomaterial for Tissue Engineering

Weizhen Sun <sup>1</sup>, David Alexander Gregory <sup>1,2</sup>, Mhd Anas Tomeh <sup>1</sup> and Xiubo Zhao <sup>1,3,\*</sup>


*Int. J. Mol. Sci.* 2021, 22, 1499. <https://doi.org/10.3390/ijms22031499>

Research Article

## Exosomes of adult human fibroblasts cultured on 3D silk fibroin nonwovens intensely stimulate neoangiogenesis


Peng Hu <sup>1,2,†</sup>, Anna Chiarini <sup>1,†,\*</sup>, Jun Wu <sup>1,3</sup>, Giuliano Freddi <sup>4</sup>, Kaiyu Nie <sup>2</sup>, Ubaldo Armato <sup>1,3</sup> and Ilaria Dal Prà <sup>1,3,†,\*</sup>

*Burns & Trauma*, Volume 9, 2021, tkab003, <https://doi.org/10.1093/burnst/tkab003>



<https://creativecommons.org/licenses/by/4.0/>

# Propagation of Neonatal Human Dermal Fibroblasts (HDFn) on Bio-compatible, Silk-derived Fibroin Discs (Mats). Cells Cultured Under Chemically-defined Animal Origin-free (cdAOF) Cell Culture Conditions for 4.5 Days, Using AvantBio's HKSdaFREE KIT Supplements, +/- Pretreatment with Recombinant Human Collagen-1




**Methods:** As described for slide4, serum-starved, post-primary passage 2 human neonatal dermal fibroblast (HDFn) were cultured under chemically defined animal origin-free (cdAOF) conditions using HFSdaFREE KIT supplements (as described in the Figure). HDFn were cultured on silk-derived Fibroin discs (mats in 24-well format), in triplicate, both with and without human recombinant Collagen-1 (rec hu Coll-1) pre-coating. Results are displayed as final mean cell densities after 4.5 days of culture, +/- the SEM. Seeding density = 5,000 cells/cm<sup>2</sup> (see red arrow). Pop. Dbl./day = average population doublings per day (a measure of relative growth rate). EpiLife Medium (ELM) and medium m106 are registered trademarks of Life Technologies Corporation.



# Propagation of Neonatal Human Dermal Fibroblasts (HDFn) on Bio-compatible, Silk-derived Fibroin Discs (Mats). Cells Cultured Under Chemically-defined Animal Origin-free (cdAOF) Cell Culture Conditions for 4.5 Days, Using AvantBio's HKSdaFREE KIT Supplements, +/- Pretreatment with Recombinant Human Collagen-1

## Photomicrographs after 1 and 4.5 days of Culture



**Methods:** As described for slide 4, serum-starved, post-primary passage 2 human neonatal dermal fibroblast (HDFn) were cultured under chemically defined animal origin-free (cdAOF) conditions using HFSdaFREE KIT supplements (as described in the previous figure). HDFn were cultured on silk-derived Fibroin discs (mats in 24-well format), in triplicate, both with and without human recombinant Collagen-1 (rec hu Coll-1) pre-coating. Results are presented as photomicrographs taken at day 1 and 4.5. EpiLife Medium (ELM) and medium m106 are registered trademarks of Life Technologies Corporation.

**WI-38 and Its Close Relative, MRC-5, are Normal Diploid Human Fetal Lung Fibroblast Lines That Have Been Utilized for Decades to Produce a Variety of Human Vaccines. WI-38 and MRC-5 Are Still Utilized Today, to Produce These Human Vaccines.**

*Over a billion vaccine units have been produced in WI-38. Examples include: poliomyelitis, measles, mumps, rubella, varicella (chicken pox), herpes zoster, adenovirus, rabies, Hepatitis A*



#### **The Role of the WI-38 Cell Strain in Saving Lives and Reducing Morbidity**


**S. J. Olshansky and L. Hayflick**

AIMS Public Health. 2017; 4(2): 127–138. Published online 2017 Mar 2.

***“Here we illustrate how the discovery and use of a single cell strain used to grow most viral vaccines in use today (WI-38 [8] and a later derivative [9]), has already had a powerful impact on human life on an order of magnitude that is unprecedented in the history of public health<sup>2</sup>.”***

***“WI-38 and its derivatives are still in use for producing many viral vaccines that are distributed worldwide today.”***

# Efficient Serial Propagation of Cryopreserved, Serum-starved WI-38 Human Diploid Fetal Lung Fibroblasts Under Chemically Defined Animal Origin-free Cell (cdAOF) Culture Conditions, Using the Experimental HFSdaFREE2 v1 KIT System Supplements




**Methods:** Cryopreserved, passage 15, population doubling level 23, WI-38 fetal human lung fibroblasts from ATCC (ATCC CCL-75), were thawed and then placed into post-primary passage 16 cell culture using EMEM + 10% FBS (from ATCC). Cells were plated into 3 T-75 flasks, at 3,000 viable cells/cm<sup>2</sup> and subsequently cultured, with 2 feedings, for 6 days. On day 6, the cells were washed 2 times over 30 minutes, with 50 mls unsupplemented EMEM basal medium (serum-starved WI-38). The WI-38 cell cultures were dissociated using chemically-defined animal origin-free (cdAOF) TrypLE Select (Life Technologies) and then cryopreserved in cdAOF CRYOVIVE 5% DMSO cryofluid. The serum-starved passage 16 WI-38 cells were later thawed and directly plated into passage 17 culture at 2,500 cells/cm<sup>2</sup> using EMEM + 10% FBS (control) or HFSdaFREE2 KIT supplements (as described in the figure) on rh Collagen-1 precoated 6 well plates, in triplicate. WI-38 were serially passaged through passage 24. Mean cumulative population doublings were calculated at the end of each passage, from triplicate wells, +/- the SEM. EpiLife and Medium (m106) are trademarks of Life Technologies Corporation.



# Efficient Serial Propagation of Cryopreserved, Serum-starved WI-38 Human Diploid Fetal Lung Fibroblasts Under Chemically Defined Animal Origin-free Cell (cdAOF) Culture Conditions, Using the Experimental HFSdaFREE2 v1 KIT System Supplements

## Photomicrographs





**Methods:** Cryopreserved, passage 15, population doubling level 23, WI-38 fetal human lung fibroblasts from ATCC (ATCC CCL-75), were thawed and then placed into post-primary passage 16 cell culture using EMEM + 10% FBS (from ATCC). Cells were plated into 3 T-75 flasks, at 3,000 viable cells/cm<sup>2</sup> and subsequently cultured, with 2 feedings, for 6 days. On day 6, the cells were washed 2 times over 30 minutes, with 50 mls unsupplemented EMEM basal medium (serum-starved WI-38). The WI-38 cell cultures were dissociated using chemically-defined animal origin-free (cdAOF) TrypLE Select (Life Technologies) and then cryopreserved in cdAOF CRYOVIVE 5% DMSO cryofluid. The serum-starved passage 16 WI-38 cells were later thawed and directly plated (without stepwise adaption to 0% FBS) into passage 17 culture at 2,500 cells/cm<sup>2</sup> using EMEM + 10% FBS (control) or HFSdaFREE2 KIT supplements (as described in the figure) on rh Collagen-1 pre-coated 6 well plates, in triplicate. WI-38 were serially passaged through passage 24. Mean cumulative population doublings (cPDL) were calculated at the end of each passage, from triplicate wells, +/- the SEM. Representative photomicrographs were taken at the day and passage level indicated in the figure. EpiLife and Medium (m106) are trademarks of Life Technologies Corporation.

# Human Hair Follicle Dermal Papilla Cells (HFDPC): Potential Use for Hair Regeneration / Restoration



# Chemically-defined Animal Origin-free (cdAOF) Propagation of Serum-Starved, Late-passage Adult\* Human Hair Follicle Dermal Papilla Cells (HFDPC), After 8 years in Cryopreservation: Comparison of AvantBio's HFSdaFREE KIT Supplement System to an Animal-originated Cell Culture (FBS) System, on Two Different Substrata (\*53 year old female temporal scalp hair, from Cell Applications, Inc.)



**Methods:** Post-primary passage 2, cryopreserved adult human dermal papilla cells (HFDPC), from Cell Applications Inc., originally established in Papilla Cell Growth Medium (FBS + growth factors), were thawed and placed directly into chemically-defined animal origin-free (cdAOF) cell culture conditions. Briefly, cryopreserved HFDPC were plated into passage 3 culture at 2,000 viable cells per cm<sup>2</sup>, under cdAOF culture conditions, into human recombinant collagen-1 precoated T-25 flasks using HFSdaFREE KIT supplements in a 50/50 blend of EpiLife medium and medium m106 (Life Technologies). HFDPC were cultured for 5 days, with feedings, dissociated with TrypLE select and then cryopreserved as passage 3 serum-starved HFDPC cells in AvantBio's cdAOF CRYOVIVE 5% DMSO cryofluid. After 8 years storage in liquid nitrogen vapor, the passage 3 serum-starved HFDPC were plated into passage 4 culture, at 1,000 viable cells/cm<sup>2</sup>, into 12- well plates, in triplicate. Cell-well substrata consisted of untreated CellBind (Corning) cell culture-ware, bovine collagen-1 pretreated conventional cell culture-ware or untreated conventional cell culture-ware, as indicated in the figure. The Passage 4 serum-starved HFDPC were cultured for 7 days. Day 7 mean cell densities and mean population doublings/day were calculated from triplicate wells, +/- the SEM.



# Chemically-defined Animal Origin-free (cdAOF) Propagation of Serum-Starved, Late-passage Adult\* Human Hair Follicle Dermal Papilla Cells (HFDPC), After 8 years in Cryopreservation: Comparison of AvantBio's HFSdaFREE KIT Supplement System to an Animal-originated Cell Culture (FBS) System, on Two Different Substrata (\*53 year old female temporal scalp hair, from Cell Applications, Inc.)

Day 7 Stained Cells



HFSdaFREE lot 108 KIT supplements on untreated CellBind Cell Culture-ware



HFSdaFREE lot 108 KIT supplements on Bovine Collagen-1 treated Cell Culture-ware



LSGS (2% FBS + Growth Factors ) on untreated Cell Culture-ware


Day 7 Photomicrographs



HFSdaFREE lot 108 KIT supplements on untreated CellBind Cell Culture-ware



HFSdaFREE lot 108 KIT supplements on Bovine Collagen-1 treated Cell Culture-ware



LSGS (2% FBS + Growth Factors ) on untreated Cell Culture-ware

**Methods:** Post-primary passage 2, cryopreserved adult human dermal papilla cells (HFDPC), from Cell Applications Inc., originally established in Papilla Cell Growth Medium (FBS + growth factors), were thawed and placed directly into chemically-defined animal origin-free (cdAOF) cell culture conditions. Briefly, cryopreserved HFDPC were plated into passage 3 culture at 2,000 viable cells per  $\text{cm}^2$ , under cdAOF culture conditions, into human recombinant collagen-1 precoated T-25 flasks using HFSdaFREE KIT supplements in a 50/50 blend of EpiLife medium and medium m106 (Life Technologies). HFDPC were cultured for 5 days, with feedings, dissociated with TrypLE select and then cryopreserved as passage 3 serum-starved HFDPC cells in AvantBio's cdAOF CRYOVIVE 5% DMSO cryofluid. After 8 years storage in liquid nitrogen vapor, the passage 3 serum-starved HFDPC were plated into passage 4 culture, at 1,000 viable cells/ $\text{cm}^2$ , into 12- well plates, in triplicate. Cell-well substrata consisted of untreated CellBind (Corning) cell culture-ware, bovine collagen-1 pretreated conventional cell culture-ware or untreated conventional cell culture-ware, as indicated in the figure. The Passage 4 serum-starved HFDPC were cultured for 7 days. Representative-stained cells as well as photomicrographs were taken at day 7. EpiLife, medium m106 and TrypLEselect are trademarks of Life Technologies Corporation. CellBind is a trademark of Corning Inc.



# Human Mesenchymal Stem Cells (HMSC)

## *Clinical Applications in Regenerative and Esthetic Medicine Using Human Mesenchymal Stem Cells (HMSC)*

- **2D Cell-based Therapies**
- **2D Cell-derived Therapeutic Products (e.g., Conditioned Medium and Exosomes)**
- **3D Reconstructed Tissue Therapies**



## *Cosmeceutical Applications Using Human Mesenchymal Stem cells*

- **Cell-derived Products (e.g., Conditioned Medium and Exosomes)**



February 2021 Frontiers in Immunology 11 DOI:10.3389/fimmu.2020.609961 Copyright © 2021 Sanchez-Castro et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY).

# A Modified HFSdaFREE - KIT Supplement System (HFSdaFREE2 v1 KIT) Efficiently Supports the Chemically-defined Animal Origin-free (cdAOF) Serial Propagation of Adult Adipose-derived Human Mesenchymal Stem Cells (HMSCs)



**Methods:** Cryopreserved post-primary passage 1 STEMPro Adipose-derived Human Mesenchymal Stem Cells (HMSCs), reared in MesenPro RS Medium (2% FBS + growth factors + basal medium) were purchased from Life Technologies. HMSCs were thawed and plated directly into passage 2 culture using MesenPro RS Medium. At the end of passage 2, the cells were washed 3 times with chemically-defined animal origin-free (cdAOF) FBM basal medium (LONZA), dissociated with the cdAOF TrypLE Select reagent (Life Technologies) and then cryopreserved using AvantBio's cdAOF CROVIVE 5% DMSO cryopreservation media. Briefly, using AvantBio's experimental cdAOF HFSdaFREE2 v1 KIT supplements (HFSdaFREE2 v1, HFGE, HFGE2), in a 50/50 blend of cdAOF KBM and FBM basal medium (LONZA), the serum-starved passage 2 HMSCs were thawed and plated into passage 3 culture at 2,500 cells/cm<sup>2</sup>. The HMSCs were serially passaged through passage 6, as shown in the figure above. Mean cumulative population doublings were calculated at the end of each passage, from triplicate wells, +/- the SEM. Representative photomicrographs at each end of passage 3 and passage 4. STEMPro cells, MesenPro RS medium and TrypLE Select are trademarks of Life Technologies Corporation. FBM and KBM basal medium are trademarks of LONZA Corporation.



AvantBio Corporation  
M: P.O. BOX 5820, Lynnwood, WA 98046  
P: [206.823.5895](tel:206.823.5895) | E: [info@avantbio.com](mailto:info@avantbio.com)  
[www.avantbio.com](http://www.avantbio.com)

Alteration of any or all of the content in this document, or incorporation of any or all of the content in this document into other publications, is strictly prohibited by AvantBio Corporation. This document is for informational use only.

Also see: <http://www.avantbio.com/legal-statement/>

AvantBio Corporation 2023, All Rights Reserved